



Anti-Sulphadoxine polyclonal antibody (DPAB-DC4787)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Acetyl-sulphadoxine 218%
Immunogen	Sulphadoxine conjugated with HSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Sulfadoxine (also spelled sulphadoxine) is an ultra-long-lasting sulfonamide previously used in combination with pyrimethamine to treat or prevent malaria. Due to high levels of resistance, its use is no longer recommended routinely. It is also used, usua
Keywords	Sulphadoxine