



Anti-Sulphamerazine polyclonal antibody (DPAB-DC4789)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | Sulphamerazine conjugated with HSA. |
| Source/Host | Sheep |
| Species Reactivity | N/A |
| Conjugate | Unconjugated |
| Applications | EIA, |
| Format | Liquid |
| Size | 100 µl |
| Buffer | In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide) |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Sulfamerazine is a sulfonamide antibacterial. |
| Keywords | Sulphamerazine |