



Anti-Sulphamerazine polyclonal antibody (DPAB-DC4790)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Acetyl Sulphamethazine 23%
Immunogen	Sulphamethazine conjugated with HSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Sulfamerazine is a sulfonamide antibacterial.
Keywords	Sulphamerazine; Sulphamethazine