



Anti-Swine influenza virus H1N1 Hemagglutinin Polyclonal antibody (DPAB-DC4798)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The peptide sequence is unique from the peptide sequence for product DPAB-DC4799 and DPAB-DC4800. This antibody is a cognate pair with product DPAB-DC4724.
Target	Swine influenza virus H1N1 Hemagglutinin
Immunogen	A synthetic peptide corresponding to Swine H1N1 Hemagglutinin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Swine influenza virus
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

The Pandemic H1N1/09 virus is a swine origin Influenza A virus subtype H1N1 virus strain responsible for the 2009 flu pandemic.

Keywords

Hemagglutinin; Swine H1N1 Hemagglutinin
