



# Anti-Swine influenza virus H1N1 Hemagglutinin Polyclonal antibody (DPAB-DC4800)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	The peptide sequence is unique from the peptide sequence for product DPAB-DC4799 and DPAB-DC4798. This antibody is a cognate pair with product number DPAB-DC4726.
<b>Target</b>	Swine influenza virus H1N1 Hemagglutinin
<b>Immunogen</b>	A synthetic peptide corresponding to Swine H1N1 Hemagglutinin.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Swine influenza virus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Introduction**

The Pandemic H1N1/09 virus is a swine origin Influenza A virus subtype H1N1 virus strain responsible for the 2009 flu pandemic.

---

**Keywords**

Hemagglutinin; Swine H1N1 Hemagglutinin

---