



Anti-Swine influenza virus Nonstructural Protein 1 Polyclonal antibody (DPAB-DC4802)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Swine influenza virus Nonstructural Protein 1
Immunogen	A synthetic peptide corresponding to Swine-Origin Influenza NS1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Swine influenza virus
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The Pandemic H1N1/09 virus is a swine origin Influenza A virus subtype H1N1 virus strain responsible for the 2009 flu pandemic.
--------------	---

Keywords

NS1; Swine-Origin Influenza NS1
