



# Anti-Testosterone (C-terminal) polyclonal antibody (DPAB-DC4810)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	4-ANDROSTEN-17b-OL-3-ONE3-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Testosterone is a steroid hormone from the androgen group and is found in mammals, reptiles, birds, and other vertebrates. In mammals, testosterone is secreted primarily by the testicles of males and the ovaries of females, although small amounts are also
Keywords	Testosterone