



Anti-Tetanus Toxoid polyclonal antibody (DPAB-DC4812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Tetanus toxoid.
Immunogen	Purified tetanus toxoid.
Source/Host	Goat
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Tetanus toxoid
----------	----------------