



Anti-TFF1 (C-terminal) polyclonal antibody (DPAB-DC4816)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a site of high antigenicity near the C-terminus of TFF1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Lyophilized
Size	500 µl
Buffer	Lyophilized from PBS (1% BSA, 15 mM sodium azide)
Preservative	15mM Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with steril distilled water, Store at 4°C for at least two months. For long term storage store at -20°C or below. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords TFF1

GENE INFORMATION

Entrez Gene ID [7031](#)

UniProt ID [P04155](#)
