



Anti-Thiabendazole polyclonal antibody (DPAB-DC4818)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Thiabendazole conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Thiabendazole (INN, BAN), thiabendazole (AAN, USAN), TBZ and the trade names Mintezol, Tresaderm, and Arbotect) is a fungicide and parasiticide.
Keywords	Thiabendazole