



Anti-Thiamine polyclonal antibody (DPAB-DC4819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Thiamine (Vitamin B1) conjugated with carrier protein.
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Thiamine or thiamin or vitamin B1, named as the "thio-vitamine" ("sulfur-containing vitamin") is a water-soluble vitamin of the B complex. First named aneurin for the detrimental neurological effects if not present in the diet, it was eventually assigned
Keywords	Vitamin B1; Thiamine (Vitamin B1)