



# Anti-Thromboxane B2 polyclonal antibody (DPAB-DC4822)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This is a sheep polyclonal antibody specific to Thromboxane B2 (TXB2). HAPTEN REACTIVITY: TXB2 100%, keto-PGF1alpha 0.17%, PGF1alpha 0.14%, PGF2alpha 0.13%, PGD2 <0.07%, 6-KetoPGF1alpha <0.07%, 13,14-dihydro-15-keto PGF2alpha <0.07%, 13,14-dihydro- 6,15-diket PG
<b>Immunogen</b>	Thromboxane B2.
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, RIA,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Thromboxane B2 is an inactive metabolite/product of thromboxane A2. It is almost completely cleared in the urine. It itself is not involved in platelet activation and aggregation in case of a wound, but its precursor, thromboxane A2, is. Thromboxane A2 sy
---------------------	---

**Keywords**Thromboxane B2

---