



# Anti-Thromboxane B2 polyclonal antibody (DPAB-DC4823)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Thromboxane B2 conjugated with KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (30% glycerol, 1% BSA, 0.02% thimerosal)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Thromboxane B2 is an inactive metabolite/product of thromboxane A2. It is almost completely cleared in the urine. It itself is not involved in platelet activation and aggregation in case of a wound, but its precursor, thromboxane A2, is. Thromboxane A2 sy
Keywords	Thromboxane B2