



# Anti-Thyroxine polyclonal antibody (DPAB-DC4825)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	T4.
Immunogen	T4/Bovine Thyroglobulin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	The thyroid hormones, triiodothyronine (T3) and its prohormone, thyroxine (T4), are tyrosine-based hormones produced by the thyroid gland that are primarily responsible for regulation of metabolism. Iodine is necessary for the production of T3 and T4. A d
--------------	---

**Keywords**

Thyroxine

---