



Anti-Tobramycin polyclonal antibody (DPAB-DC4828)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Tobramycin conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Tobramycin is an aminoglycoside antibiotic derived from <i>Streptomyces tenebrarius</i> and used to treat various types of bacterial infections, particularly Gram-negative infections. It is especially effective against species of <i>Pseudomonas</i> .
Keywords	Tobramycin