



Anti-T. gondii Polyclonal antibody (DPAB-DC4830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	T. gondii
Immunogen	Toxoplasma gondii
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. gondii
Purification	Immunoaffinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	In 10 mM PBS, pH 7.4 (1% BSA, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Toxoplasma gondii is an obligate, intracellular, parasitic protozoan that causes the disease toxoplasmosis.
Keywords	Toxoplasma Gondii