



# Rabbit Anti-T. gondii trophozoites Polyclonal antibody (DPAB-DC4831)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Toxoplasma gondii trophozoites. Cross-reactivity to tachyzoites has not been tested.
<b>Target</b>	T. gondii
<b>Immunogen</b>	RH strain of T. gondii, purified trophozoites.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	T. gondii
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Toxoplasma gondii is an obligate, intracellular, parasitic protozoan that causes the disease
---------------------	--

toxoplasmosis.

---

**Keywords**

Toxoplasma Gondii

---