



Anti-T. gondii TGME49_079100 Polyclonal antibody (DPAB-DC4832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	T. gondii TGME49_079100
Immunogen	A synthetic peptide corresponding to internal region of Toxoplasma gondii protein TGME49_079100.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. gondii
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Toxoplasma gondii is an obligate, intracellular, parasitic protozoan that causes the disease toxoplasmosis.
--------------	---

Keywords

protein TGME49_079100; Toxoplasma gondii protein TGME49_079100
