



Anti-Trans-hydroxyproline polyclonal antibody (DPAB-DC4835)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity ratioCis-Hydroxyproline-G-carrier protein : 1
Immunogen	Trans-hydroxyproline conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere.After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years.Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Trans-hydroxyproline