



Anti-Trenbolone (C-terminal) polyclonal antibody (DPAB-DC4836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	4, 9, 11-ESTRATRIEN-17b-0L-3-ONE3-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate
Keywords	Trenbolone