



Anti-T. pallidum Polyclonal antibody (DPAB-DC4837)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes all antigens. Non-reactive in RPR assay. Has not been absorbed and may react with related microorganisms.
Target	T. pallidum
Immunogen	Native Treponema pallidum.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. pallidum
Conjugate	Unconjugated
Applications	IHC-P, IF,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Treponema pallidum is a spirochaete bacterium with subspecies that cause treponemal diseases such as syphilis, bejel, pinta, and yaws. The treponemes have a cytoplasmic and an outer membrane. Using light microscopy, treponemes are only visible using dark

Keywords

Treponema pallidum
