



Anti-Triclabendazole polyclonal antibody (DPAB-DC4838)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Triclabendazole conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Triclabendazole (commercial names: veterinary, liquid: Fasinex; human, tablets: Egaten, both manufactured by Novartis) is a member of the benzimidazole family of anthelmintics. The benzimidazole drugs share a common molecular structure, triclabendazole be
Keywords	Triclabendazole