



Anti-Troponin I-Cardiac polyclonal antibody (DPAB-DC4840)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific for human Cardiac Troponin I. Recognizes epitopes located in the 27-39 region of Troponin I. Cross reactivity with human skeletal Troponin I is 0.02% by ELISA.
Immunogen	Purified human cardiac Troponin I.
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC, ELISA, sELISA, Detection Antibody,
Format	Liquid
Size	1 mg
Buffer	In 10 mM sodium phosphate buffer, pH 7.2 (0.85% NaCl, 0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Keywords Troponin I-Cardiac

GENE INFORMATION

Entrez Gene ID [7139](#)

UniProt ID [P45379](#)
