



Anti-TSST-1 polyclonal antibody (DPAB-DC4843)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Toxic Shock Syndrome Toxin. No reactivity with other Staph toxins.
Immunogen	Purified TSST-1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ID,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	TSST-1
----------	--------