



## Anti-Tub-1 (aa 160-421) polyclonal antibody (DPAB-DC4845)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 160-421 of tub-1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Nematoda
Conjugate	Unconjugated
Applications	WB, IP,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C or -80°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

Keywords	Tub-1
----------	-------