



# Anti-Tubulin-alpha (C-terminal) polyclonal antibody (DPAB-DC4846)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of human Tubulin-alpha.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In HEPES, 150 mM NaCl (50% glycerol, 0.01% BSA, 0.03% sodium azide)
<b>Preservative</b>	0.03% Sodium Azide
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	Tubulin-alpha
-----------------	---------------