



Anti-Ubiquitin (native) polyclonal antibody (DPAB-DC4850)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified bovine Ubiquitin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Chicken, Dog, Frog, Guinea pig, Hamster, Human, Monkey, Mouse, Pig, Rabbit, Rat, Sheep, Yeast
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	50 µg
Buffer	In PBS
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ubiquitin is a small (8.5 kDa) regulatory protein that has been found in almost all tissues (ubiquitously) of eukaryotic organisms. It was discovered in 1975 by Goldstein and further characterized throughout the 1970s and 1980s. There are four genes in th
Keywords	Ubiquitin