



Anti-Ubiquitin polyclonal antibody (DPAB-DC4853)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	human ubiquitin protein.
Immunogen	Purified ubiquitin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-F, ICC/IF, IP
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.6, with 1% BSA
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ubiquitin is a small (8.5 kDa) regulatory protein that has been found in almost all tissues (ubiquitously) of eukaryotic organisms. It was discovered in 1975 by Goldstein and further characterized throughout the 1970s and 1980s. There are four genes in th
--------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Keywords

Ubiquitin
