



Anti-Vaccinia virus G9R Polyclonal antibody (DPAB-DC4861)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Vaccinia virus G9R
Immunogen	A synthetic peptide corresponding to N-terminus of Vaccinia virus G9R.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Vaccinia virus
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Vaccinia virus (VACV or VV) is a large, complex, enveloped virus belonging to the poxvirus family. It has a linear, double-stranded DNA genome approximately 190 kbp in length, and which encodes approximately 250 genes. The dimensions of the virion are rou
--------------	---

Keywords

G9R; Vaccinia virus G9R
