



# Anti-Ariola virus Variola virus Protein J5 Polyclonal antibody (DPAB-DC4868)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Ariola virus Variola virus Protein J5
Immunogen	A synthetic peptide corresponding to C-terminus of Variola virus Protein J5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Ariola virus Variola virus
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Smallpox is caused by the variola virus that emerged in human populations thousands of years ago. Except for laboratory stockpiles, the variola virus has been eliminated. However, in the aftermath of the events of September and October, 2001, there is hei
--------------	---

**Keywords**

Protein J5; Variola virus Protein J5

---