



Anti-Virginiamycin M polyclonal antibody (DPAB-DC4874)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Virginiamycin M conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Virginiamycin is a streptogramin antibiotic similar to pristinamycin and quinupristin/dalfopristin. It is a combination of pristinamycin IIA (virginiamycin M1) and virginiamycin S1. Virginiamycin is used in the fuel ethanol industry to prevent microbial c
Keywords	Virginiamycin M