



# Anti-WNV Core Protein Polyclonal antibody (DPAB-DC4877)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	WNV Core Protein
<b>Immunogen</b>	A synthetic peptide corresponding to C-terminus 15 amino acids of West Nile Virus Core protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	WNV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Transfected lysate), ICC, IF, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	West Nile virus (WNV) is a mosquito-borne zoonotic arbovirus belonging to the genus Flavivirus in the family Flaviviridae. This flavivirus is found in temperate and tropical regions of the world.
---------------------	---

It was first identified in the West Nile subregion in the E

---

**Keywords**

Core; West Nile Virus Core

---