



Anti-WNV Envelope Antigen Polyclonal antibody (DPAB-DC4878)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	WNV Envelope Antigen
Immunogen	A synthetic peptide corresponding to internal region 15 amino acids of West Nile virus Envelope.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	WNV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	West Nile virus (WNV) is a mosquito-borne zoonotic arbovirus belonging to the genus Flavivirus in the family Flaviviridae. This flavivirus is found in temperate and tropical regions of the world.
---------------------	---

It was first identified in the West Nile subregion in the E

Keywords

Envelope; West Nile Virus Envelope
