



Anti-Y. pseudotuberculosis Polyclonal antibody (DPAB-DC4883)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Y. pseudotuberculosis
Immunogen	Yersinia pseudotuberculosis
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Y. pseudotuberculosis
Conjugate	Unconjugated
Applications	ICC, IF,
Format	Lyophilized
Size	200 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with 200 uL deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Yersinia pseudotuberculosis is a Gram-negative bacterium that causes Far East scarlet-like fever in humans, who occasionally get infected zoonotically, most often through the food-borne route. Animals are also infected by Y. pseudotuberculosis. The bacter
--------------	---

Keywords

Y. pseudotuberculosis; Yersinia pseudotuberculosis
