



Anti-C21ORF128 polyclonal antibody (DPAB-DC663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C21orf128 (chromosome 21 open reading frame 128) is a protein-coding gene, and is affiliated with the lncRNA class.
Immunogen	Recombinant protein corresponding to amino acids of recombinant C21orf128. The sequence is GLPCHQNTAGANPHLFLGCTSSSLQGLEYGGQRGDAHGKPGVLHGELEPHDHTSRLERH DLHSQLPTSVQVRHHWWE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C21orf128 chromosome 21 open reading frame 128 [Homo sapiens (human)]
Official Symbol	C21ORF128
Synonyms	C21ORF128; chromosome 21 open reading frame 128;
Entrez Gene ID	150147
Chromosome Location	21q22.3