



Anti-HC polyclonal antibody (DPAB-DC669)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recombinant mouse C5a can be used as a positive control.
Immunogen	Recombinant protein corresponding to mouse Hc.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS pH 7.2 (30% glycerol, 0.5% BSA, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Hc hemolytic complement [Mus musculus (house mouse)]
Official Symbol	HC
Synonyms	HC; hemolytic complement; C5; He; C5a; complement C5;

Entrez Gene ID	15139
Protein Refseq	NP_034536
UniProt ID	P06684
Chromosome Location	2 cen-C1; 2 23.22 cM
Pathway	Activation of C3 and C5; Complement Activation, Classical Pathway; Complement and coagulation cascades; Complement cascade.
Function	endopeptidase inhibitor activity;