



Anti-DTX3L (aa 3-110) polyclonal antibody (DPAB-DC672)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DTX3L functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).[supplied by OMIM, Nov 2009]
Immunogen	DTX3L (NP_612144, 3 a.a. ~ 110 a.a) partial recombinant protein with GST tag. The sequence is SHLRPPSPLLVRVYKSGPRVRRKLESYFQSSKSSGGGECTVSTQEHEAPGTRVEFSERA AKERVLKKGEHQILVDEKPVPIFLVPTENSIKKNTRPQISSLTQSQAE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DTX3L deltex 3 like, E3 ubiquitin ligase [Homo sapiens (human)]
Official Symbol	DTX3L

Synonyms	DTX3L; deltex 3 like, E3 ubiquitin ligase; BBAP; E3 ubiquitin-protein ligase DTX3L; rhysin2; rhysin 2; rhysin-2; deltex 3-like; B-lymphoma- and BAL-associated protein;
Entrez Gene ID	151636
Protein Refseq	NP_612144
UniProt ID	Q8TDB6
Chromosome Location	3q21.1
Pathway	Adaptive Immune System; Proteasome degradation; presentation; Notch Signaling Pathway
Function	histone binding; ligase activity; protein binding; ubiquitin-protein transferase activity