



Anti-SH3RF2 (aa 640-729) polyclonal antibody (DPAB-DC678)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SH3RF2 (SH3 domain containing ring finger 2) is a protein-coding gene. GO annotations related to this gene include protein phosphatase inhibitor activity and ligase activity. An important paralog of this gene is SH3RF3.
Immunogen	SH3RF2 (NP_660205, 640 a.a. ~ 729 a.a) partial recombinant protein with GST tag. The sequence is VKT VRFQNYSPPTKHYSHTSGKPEQPATLKASQPEAASLGPEMTVLFAHRSGCHSGQ QTDLRK SALAKATTLVSTASGTQTVFPSK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [SH3RF2 SH3 domain containing ring finger 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	SH3RF2
Synonyms	SH3RF2; SH3 domain containing ring finger 2; HEPP1; POSHER; RNF158; PPP1R39; putative E3 ubiquitin-protein ligase SH3RF2; RING finger protein 158; POSH-eliminating RING protein; SH3 domain-containing RING finger protein 2; heart protein phosphatase 1-binding protein; protein phosphatase 1 regulatory subunit 39; protein phosphatase 1, regulatory subunit 39;
Entrez Gene ID	153769
Protein Refseq	NP_689763
UniProt ID	Q08AM8
Chromosome Location	5q32
Function	ligase activity; phosphatase binding; protein binding; protein phosphatase 1 binding