



Anti-RNF145 (aa 592-691) polyclonal antibody (DPAB-DC679)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF145 (ring finger protein 145) is a protein-coding gene. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is SYVN1.
Immunogen	FLJ31951 (NP_653327, 592 a.a. ~ 691 a.a) partial recombinant protein with GST tag. The sequence is WLYVQETCP LCHCHLKNSSQLPGLGTEPVLQPHAGAEQNVMFQEGTEPPGQEHTPGTRIQ EGSRDNNEYIARRPDNQEGAFDPKEYPHSAKDEAHPVESA
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNF145 ring finger protein 145 [Homo sapiens (human)]
Official Symbol	RNF145

Synonyms	RNF145; ring finger protein 145; RING finger protein 145;
Entrez Gene ID	153830
Protein Refseq	NP_001186309
UniProt ID	Q8NDT8
Chromosome Location	5q33.3
Function	zinc ion binding;