



Anti-FMR1NB (full length) polyclonal antibody (DPAB-DC691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FMR1NB (fragile X mental retardation 1 neighbor) is a protein-coding gene. Diseases associated with FMR1NB include dyslexia, and sarcoma.
Immunogen	FMR1NB (AAH34320, 1 a.a. ~ 255 a.a) full-length recombinant protein with GST tag. The sequence is MSSHRRKAKGRNRRSHRAMRVAHLELATYELAATESNPESHHPGYEAAMADRPQPGWRES LKMRVSKPFGMLMLSIWILLFVCYYLSYYLCSGSSYFVLANGHILPNSENAHGQSLEEDS ALEALLNFFFFPTTCNLRENQVAKPCNELQDLSESECLRHKCCFSSSGTTSFKCFAPFRDV PKQMMQMFGGLGAISLILVCLPIYCRSLFWRSEPADDLQRQDNRVVTGLKKQRRKRKRKSE MLQKAARGREEHGDE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FMR1NB fragile X mental retardation 1 neighbor [Homo sapiens (human)]
Official Symbol	FMR1NB
Synonyms	FMR1NB; fragile X mental retardation 1 neighbor; CT37; NYSAR35; NY-SAR-35; fragile X mental retardation 1 neighbor protein; cancer/testis antigen 37; sarcoma antigen NY-SAR-35;
Entrez Gene ID	158521
Protein Refseq	NP_689791
UniProt ID	Q8N0W7
Chromosome Location	Xq28