



Anti-NIM1K (aa 342-435) polyclonal antibody (DPAB-DC730)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MGC42105 (NP_699192, 342 a.a. ~ 435 a.a) partial recombinant protein with GST tag. The sequence is KHLSETSTLKEEENEVKSTLEHLGITEEHIRNNQGRDARSSITGVYRIILHRVQRKKALE SVPVMMLPDPKERDLKKGSRVYRGIRHTSKFCSI
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NIM1K NIM1 serine/threonine protein kinase [Homo sapiens (human)]
Official Symbol	NIM1K
Synonyms	NIM1K; NIM1 serine/threonine protein kinase; NIM1; serine/threonine-protein kinase NIM1; NIM1 serine/threonine-protein kinase;

Entrez Gene ID	167359
Protein Refseq	NP_699192
UniProt ID	Q8IY84
Chromosome Location	5p12
Function	ATP binding; magnesium ion binding; protein serine/threonine kinase activity;