



Anti-CNBD1 polyclonal antibody (DPAB-DC741)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CNBD1 (cyclic nucleotide binding domain containing 1) is a protein-coding gene.
Immunogen	Recombinant protein corresponding to amino acids of recombinant CNBD1. The sequence is AEMYLPYSYDSMLSKWSTFGTLEVMPQNESETQMFSVVTEDDCEILKIPAKGYAKIKEEKIKLENMQKLKLRMCPPYYEEWPTLSIYELIALLKW
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CNBD1 cyclic nucleotide binding domain containing 1 [Homo sapiens (human)]
Official Symbol	CNBD1
Synonyms	CNBD1; cyclic nucleotide binding domain containing 1; cyclic nucleotide-binding domain-containing protein 1;
Entrez Gene ID	168975
Protein Refseq	NP_775809
UniProt ID	Q8NA66
Chromosome Location	8q21.3