



Anti-SNX31 (aa 242-341) polyclonal antibody (DPAB-DC743)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SNX31 (sorting nexin 31) is a protein-coding gene. Diseases associated with SNX31 include systemic lupus erythematosus, and lupus erythematosus. GO annotations related to this gene include phosphatidylinositol binding. An important paralog of this gene is SNX27.
Immunogen	MGC39715 (AAH31260, 242 a.a. ~ 341 a.a) partial recombinant protein with GST tag. The sequence is QNLELRFQYSEDSCWQWFVIYTKQAFLLSSCLKKMISEKMKVLAAGENTEMQIEVPEQSKS KKYHIQQSQKDYSSFLSRKSKIKIAKGDCVFGNIKEEDL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [SNX31 sorting nexin 31 \[Homo sapiens \(human\) \]](#)

Official Symbol	SNX31
Synonyms	SNX31; sorting nexin 31; sorting nexin-31;
Entrez Gene ID	169166
Protein Refseq	NP_689841
UniProt ID	Q8N9S9
Chromosome Location	8q22.3
Function	phosphatidylinositol binding;