



Anti-ATXN10 (C-terminal) polyclonal antibody (DPAB-DC749)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	human homolog is a region of a chromosome highly linked to autosomal dominant cerebellar ataxias [RGD, Feb 2006]
Immunogen	A synthetic peptide corresponding to C-terminus of rat Atxn10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Atxn10 ataxin 10 [Rattus norvegicus (Norway rat)]
Official Symbol	ATXN10

Synonyms	ATXN10; ataxin 10; Sca10; ataxin-10; ataxin 10 homolog; spinocerebellar ataxia 10 homolog; neuronal beta-catenin-like protein; spinocerebellar ataxia type 10 protein homolog;
Entrez Gene ID	170821
Protein Refseq	NP_579847
UniProt ID	Q9ER24
Chromosome Location	7q34
Function	identical protein binding; protein binding;