



Anti-SULF1 polyclonal antibody (DPAB-DC754)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	displays arylsulfatase and/or hydrolase enzyme activity; may catalyze the desulfation of glycosaminoglycans [RGD, Feb 2006]
Immunogen	Recombinant protein corresponding to rat Sulf1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB (Transfected lysate),
Format	Liquid
Size	50 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Sulf1 sulfatase 1 [Rattus norvegicus (Norway rat)]
Official Symbol	SULF1
Synonyms	SULF1; sulfatase 1; extracellular sulfatase Sulf-1; RSulfFP1; sulfatase FP;

Entrez Gene ID	171396
Protein Refseq	NP_599205
UniProt ID	Q8VI60
Chromosome Location	5q11
Function	N-acetylglucosamine-6-sulfatase activity; alkyl sulfatase activity; arylsulfatase activity; calcium ion binding