



Anti-MCPT9 polyclonal antibody (DPAB-DC759)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to mouse Mcpt9/Cma2. The sequence is CKIFKHYKDSLQ
Source/Host	Goat
Species Reactivity	Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In Tris saline, pH 7.3 (0.5% BSA, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Mcpt9 mast cell protease 9 [Mus musculus (house mouse)]
Official Symbol	MCPT9
Synonyms	MCPT9; mast cell protease 9; MMCP9; MMCP-1; MMCP-9; mast cell-specific serine protease;

Entrez Gene ID	17232
Protein Refseq	NP_034912
UniProt ID	Q35164
Chromosome Location	14 C3; 14 28.19 cM
Function	catalytic activity; hydrolase activity; peptidase activity; serine-type endopeptidase activity