



Anti-NIN polyclonal antibody (DPAB-DC780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NIN (ninein (GSK3B interacting protein)) is a protein-coding gene. Diseases associated with NIN include geniculate herpes zoster, and herpes zoster. GO annotations related to this gene include GTP binding and calcium ion binding. An important paralog of this gene is NINL.
Specificity	Specific to recombinant protein GX1316. This antibody detects mNIN protein. It also recognizes human NIN protein.
Immunogen	Recombinant GST fusion protein corresponding to 242 amino acids of mouse Nin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [Nin ninein \[Mus musculus \(house mouse\) \]](#)

Official Symbol	NIN
Synonyms	NIN; ninein; AI385615; AU024711; mKIAA1565; 3110068G20Rik;
Entrez Gene ID	18080
Protein Refseq	NP_001074922
UniProt ID	E9Q4S3
Chromosome Location	12 C3; 12
Function	GTP binding; calcium ion binding; kinase binding; nucleotide binding