



# Anti-PLSCR2 polyclonal antibody (DPAB-DC805)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	PLSCR2 (phospholipid scramblase 2) is a protein-coding gene. Diseases associated with PLSCR2 include endotheliitis. GO annotations related to this gene include phospholipid scramblase activity and calcium ion binding. An important paralog of this gene is PLSCR1.
<b>Immunogen</b>	A synthetic peptide corresponding to mouse Plscr2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Cell lysate),
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.1% proclin, 2.0% Block Ace)
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">Plscr2 phospholipid scramblase 2 [ Mus musculus (house mouse) ]</a>
------------------	---

<b>Official Symbol</b>	PLSCR2
<b>Synonyms</b>	PLSCR2; phospholipid scramblase 2; PL scramblase 2; ca(2+)-dependent phospholipid scramblase 2;
<b>Entrez Gene ID</b>	<a href="#">18828</a>
<b>Protein Refseq</b>	<a href="#">NP_001182013</a>
<b>UniProt ID</b>	<a href="#">Q9DCW2</a>
<b>Chromosome Location</b>	9 E3.3; 9
<b>Function</b>	SH3 domain binding; metal ion binding; protein binding;