



Anti-PODXL polyclonal antibody (DPAB-DC816)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	inhibits cell-cell adhesion; may play a role in maintaining filtration in the glomerular epithelium [RGD, Feb 2006]
Immunogen	A synthetic peptide corresponding to rat Podxl.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB (Tissue lysate),
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Podxl podocalyxin-like [Rattus norvegicus (Norway rat)]
Official Symbol	PODXL

Synonyms	PODXL; podocalyxin-like; PC; PCLP-1; podocalyxin; podocalyxin; podocalyxin-like protein 1;
Entrez Gene ID	192181
Protein Refseq	NP_620203
UniProt ID	Q9WTQ2
Chromosome Location	4q22
Function	protein binding;