



Anti-DTX3 (internal region) polyclonal antibody (DPAB-DC832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DTX3 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).[supplied by OMIM, Nov 2009]
Immunogen	A synthetic peptide corresponding to amino acids at internal region of human DTX3. The sequence is C-KNKVTVSKPVWD
Source/Host	Goat
Species Reactivity	Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In 0.5 mg/mL Tris saline, pH 7.3 (0.02% sodium azide, 0.5% BSA)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DTX3 deltex 3, E3 ubiquitin ligase [Homo sapiens (human)]
Official Symbol	DTX3
Synonyms	DTX3; deltex 3, E3 ubiquitin ligase; RNF154; deltex3; probable E3 ubiquitin-protein ligase DTX3; deltex homolog 3; protein deltex-3; RING finger protein 154;
Entrez Gene ID	196403
Protein Refseq	NP_001273174
UniProt ID	Q8N9I9
Chromosome Location	12q13.3
Pathway	Notch Signaling Pathway;
Function	ligase activity; zinc ion binding;